

TWIST+ LICK™ ORAL GEL

BRUSH FREE™



features & benefits



ProLong Technology™
Continuous release



Cleans,
whitens &
brightens



Helps control
plaque & tartar
buildup



No brushing
(for real!)



Kissable
for days

**TWIST
& LICK**



TWIST.

Turn dial until gel appears, then click the appropriate dose.




LICK.

Dogs will eagerly lick the gel off the applicator.



DONE.

Gel bonds to gums & creates protective biofilm on gums. ProLong Technology slowly releases active ingredients helping restore good oral hygiene.



Small Dogs
< 25 lbs
2 clicks
 up to 2-month supply

Large Dogs
25 lbs & over
4 clicks
 up to 1-month supply

INGREDIENTS: Chlorhexidine, Sodium Bicarbonate, ProLONG Technology™ (proprietary blend of Sorbitol, Glycerin, Water, Soy Lecithin, Mineral Oil, Natural Chicken Flavor)